


## IMG UI Preferences

IMG UI sets a limitation on numbers of genome, genes, etc. to be display in a table display in order to improve the UI performance. Such limitations are needed in order not to overwhelm the web browser. Default setting also hides viruses, plasmids and genome fragments in the UI display. The default settings should work well for most users. However, users can overwrite the default setting using the **Preferences** option under **MyIMG** menu.

### Preferences

Parameter	Current Setting
Max. Paralog Groups	500 ▾
Max. Gene / Scaffold List Results	1000 ▾
Max. Homolog Results	200 ▾
Max. Taxon Gene Neighborhoods	5 ▾
Min. Homolog Percent Identity	30 ▾
Hide Viruses From Genome Lists	Yes ▾
Hide Plasmids From Genome Lists	Yes ▾
Hide GFragment From Genome Lists	Yes ▾
Hide Zeroes in Genome Statistics	Yes ▾
Session Cache On (Yes recommended)	Yes ▾

Save Preferences
Default Settings


**hint:** For faster processing, adjust to lower numbers. For more complete result lists, adjust to higher numbers. The default settings should work well for most users.

**Figure 1.** IMG Preferences

To change a setting, simply select a different value in the **Current Setting** drop down list and click the **Save Preferences** button. Clicking the **Default Settings** button will return to the default preference setting.